

# Phys 4910 Spectroscopy

## Monochromator Resolution

### Introduction

In the previous lab you looked at the D lines of Sodium and Potassium. For both the elements the spin-orbit coupling is reasonably large, and it is easy to see the D as separate. In the case of lithium the splitting is much smaller, and for hydrogen it is too small for us to see.

In order to be able to resolve the D lines for potassium the resolving power of the monochromator should be as high as possible. This lab has two parts; first to measure the resolving power as a function of the width of the entrance and exit slits, and secondly to try to resolve the D lines of potassium and measure the corresponding spin-orbit coupling.

The portion that you need to record is just that wide enough to include both the sodium D lines<sup>(1)</sup>. They are approximately  $17 \text{ cm}^{-1}$  apart<sup>(2)</sup>, which gives you the scale for your spectrum. You don't need the neon lines for calibration.

The entrance and exit slit widths are controlled by the micrometer located directly above the appropriate slit. Each marking on the micrometer corresponds to  $10 \mu$  of slit width. One full revolution of the micrometer therefore corresponds to changing the slit width by  $250 \mu$ .

### Definitions

- Resolving power:  $\lambda/\delta\lambda = \nu/\delta\nu$  where  $\delta\lambda$  is the minimum separation of two lines which can be distinguished as separate. We shall set this equal to the FWHM of one peak.
- Full width half maximum (FWHM): the width (measured as  $\delta\lambda$  or as  $\delta\nu$ ) at an intensity equal to  $1/2$  of the peak intensity
- Dispersion: the wavelength range (usually in nm) per unit of measurement.

### Assignment

#### *Instrument Resolving Power*

1. Position the neon/sodium lamp in front of the entrance slit, and adjust its position to get a good signal. (A neon line would work for this part.)
2. You will be repeatably recording the same portion of the spectrum, but with different sizes for the two slits. At each setting keep the width of the exit slit equal to that of the entrance slit.
3. Check that the reading on the micrometer is 1, that is the slit width is  $10 \mu$ .
4. Record a short section of the spectrum, a portion just wide enough to include both the lines that

---

1 Although the D lines have wavelengths of approximately 589.0 and 589.6 nm, the calibration of the instrument is currently about 1 nm off. Expect to see the lines at displayed wavelengths of 588.0 and 588.6 nm.

2 A more accurate number can be obtained from the energy levels of Na I, as the difference in energy of the first two excited levels, see [http://physics.nist.gov/PhysRefData/ASD/levels\\_form.html](http://physics.nist.gov/PhysRefData/ASD/levels_form.html)

you selected. Record the spectrum slowly (say 5 minutes) so that the speed of the recording equipment does not affect your results. Also, sample the data as fast as possible. Since you are only recording a short section of the whole spectrum you can afford to sample quite quickly, ten times a second, even twenty times a second.


5. Repeat step 4 for slit widths between  $10\ \mu$  (micrometer reading = 1) to  $500\ \mu$  (micrometer reading = 50, corresponding to two full turns of the micrometer).
  - a. Suggested slit widths,  $10\ \mu$ ,  $20\ \mu$ ,  $50\ \mu$ ,  $100\ \mu$ ,  $150\ \mu$ ,  $200\ \mu$ ,  $250\ \mu$ ,  $350\ \mu$ , and  $500\ \mu$ .
  - b. As you increase the entrance slit width you will let in more light, and the peaks in your spectrum will get more intense. You might have to change either the voltage applied to the photomultiplier, or the sensitivity of the electrometer, or both.

### *Spectrum of Lithium*

6. Replace the neon/sodium lamp with a neon/ lithium lamp. Like sodium, lithium has a pair of D lines, but with shorter wavelength, around  $6707\ \text{\AA}$ . Their separation is quite small. Try to resolve them as separate lines with slit widths of  $10\ \mu$ ,  $20\ \mu$ , and  $50\ \mu$ . (There are also pairs of lines at approximately  $4602\ \text{\AA}$  and  $6103\ \text{\AA}$ , which you might also try to resolve as separate lines. Note that each of these pairs is usually significantly less intense than the pair at  $6707\ \text{\AA}$ .)
7. If you scan from roughly  $665\ \text{nm}$  to  $675\ \text{nm}$  you should also get two strong neon lines for calibration. Widen the scan a little if you don't get both Ne lines.

### Sample Spectrum

Please refer to this diagram for a definition of terms.



### Analysis

You will need to make a chart for each of the spectra (similar to the diagram above), one for each slit width and one per page. Then for each spectrum

1. Measure the distance between peaks (center to center). In sodium this corresponds to a spectral difference of  $17 \text{ cm}^{-1}$ .
2. Measure the width of the larger peak at the position half way down (the FWHM).  
Note: if your spectrum has a significant amount of noise then this is not at the intensity which is half the peak intensity. For example, if your peak has a peak intensity of 7000, but there is a noise (background) value of 1000, then the peak is 6000 high, and half of this is 3000. You should therefore measure the width at an intensity value of 4000, not at 3500 ( $\frac{1}{2}$  the peak value.)
3. Knowing the separation of the peaks, their difference in wavelength, and the width of the taller peak, calculate its width in terms of wavelength. Use that to obtain the resolution.
4. When you have analyzed all the spectra, plot the resolution as a function of slit width.

## Report

For inclusion in this week's lab I would like the following:

1. Abstract
2. Introduction. Most of this would be a discussion of the variation of the strength of spin orbit coupling with atomic number.
3. A brief description of your experiment.
4. One (but only one) representative graph of your Na spectrum. (Pick the best looking one!) There is no need to give me the raw data.
5. The plot of FWHM vs entrance slit width. (Again, I don't need to see raw data.)
6. Your result for the spin orbit coupling in lithium
7. A discussion of results.
8. Conclusions